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References (from Editorial)

1. Effectiveness of introducing blood culture collection packs to reduce contamination rates, A.I. Bamber, J.G. Cunniffe, D. Nayar, R. Ganguly, E. Falconer, *British Journal of Biomedical Science* 2009 66 (1)
2. Controlling blood-culture contamination rates, Dennis J. Ernst (ASCP), MLO, March 2004
3. Investigation of Blood Cultures (for organisms other than mycobacterium species), Issue n°. 6, dd. 06.07.10, Issued by: Standards Unit, Department for Evaluations, Standards and Training, Page n°. 1 of 34, BSOP37i6
4. Taking blood cultures, a summary of best practice, NHS, Produced by COI for the Department of Health, 2007
5. Cumitech 1C: Blood Cultures IV, Baron et All, ASM Press, 2005
6. Effects of Volume and Periodicity on Blood Cultures, James Li, James J. Plorde, Larry G. Carlson, *Journal of Clinical Microbiology*, Nov 1994, p. 2829-2831
7. Optimal Testing Parameters for Blood Cultures, F. R. Cockerill III, J. W. Wilson, E. A. Vetter, K. M. Goodman, C. A. Torgerson, W. S. Harmsen, C. D. Schleck, D. M. Ilstrup, J. A. Washington II, and W. R. Wilson, *Clinical Infectious Diseases*
8. Detection of Bloodstream Infections in Adults: How Many Blood Cultures Are Needed? Andrew Lee, Stanley Mirrett, L. Barth Reller, Melvin P. Weinstein, *Journal of Clinical Microbiology*, Nov. 2007 p. 2546-3548
9. Is the Volume of Blood Cultured Still a Significant Factor in the Diagnosis of Bloodstream Infections? Emilio Bouza, Dolores Sousa, Marta Rodríguez-Créixems, Juan García Lechuz, Patricia Muñoz *JOURNAL OF CLINICAL MICROBIOLOGY*, Sept. 2007, p. 2765–2769
10. Effects of Volume and Site of Blood Draw on Blood Culture Results, Wilson I Gonzales, Nancy Cornish, Michael Morre, Aimin Chen and Meera Varman, *J. Clin. Microbiol.* 2009, 47(11):3482. Sept. 2009
11. Frequency of Low-Level Bacteremia in Children from Birth to Fifteen Years of Age, JAMES A. KELLOGG, JOHN P. MANZELLA, DAVID A. BANKERT, *J. Clin. Microbiol.* June 2000, p. 2181-2185
12. How Reliable Is a Negative Blood Culture Result? Volume of Blood Submitted for Culture in Routine Practice in a Children's Hospital, Thomas G. Connell, Mhisti Rele, Donna Cowley, Jim P. Buttery and Nigel Curtis, *PEDIATRICS*, Official Journal of the American Academy of Pediatrics, 2007; 119; 891-896
13. Package Inser BacT/ALERT® SA, bioMérieux, (c) BIOMÉRIEUX 2000

VACUETTE® Safety Blood Collection Set + Blood Culture Holder

The Safety Blood Collection Set with pre-attached holder has been designed for use with most commonly used Blood Culture Bottles and **VACUETTE®** blood collection tubes in one procedure. The set is sterile and pre-assembled for ease of use, cost effective and ensures your patient only has to undergo one venipuncture for all procedures.

In order to provide a complete product line for all requirements, a standard Blood Collection Set (without safety mechanism) is also available with pre-assembled Blood Culture Holder.

Suitable for BacT / Alert and BacTec blood culture bottles.

Instructions for Use

Version a = Safety Blood Collection Set; Version b = Blood Collection Set

1 Prepare venipuncture site and clean septum of Blood Culture Bottles in accordance with your institution's protocol

2 Using the label scale as your guide, mark on the Blood Culture Bottle the required fill level


3a Perform venipuncture with **VACUETTE®** Safety Blood Collection Set

3b Perform venipuncture with **VACUETTE®** Blood Collection Set

4 Push the holder down over the aerobic Blood Culture Bottle, ensuring holder and bottle are below the puncture site. Once filled to required level, remove and repeat with anaerobic bottle. Continue with **VACUETTE®** Blood Collection Tubes (if required)

5a After completion of blood collection, hold the gauze down with one hand. With the other hand, activate the safety mechanism by pressing in both sides of the hub to engage the lock. Slide the safety mechanism backward until an audible click is heard. The click is a sign that the safety mechanism has been correctly activated

5b After completion of venipuncture, gently remove needle from vein, applying pressure to puncture site with dry sterile swab until bleeding stops

6 Promptly dispose of in an approved sharps disposal container in accordance with the procedures of your facility

Item Numbers

Item No.	Description	Size	Packaging	
			Inner	Outer
VACUETTE® Safety Blood Collection Sets + Blood Culture Holder (sterile)				
450182	Safety Blood Collection Set + Blood Culture Holder	tubing length: 19cm 21G x 3/4"	24 pcs.	240 pcs.
450183	Safety Blood Collection Set + Blood Culture Holder	tubing length: 19cm 23G x 3/4"	24 pcs.	240 pcs.
450184	Safety Blood Collection Set + Blood Culture Holder	tubing length: 30cm 21G x 3/4"	24 pcs.	240 pcs.
450185	Safety Blood Collection Set + Blood Culture Holder	tubing length: 30cm 23G x 3/4"	24 pcs.	240 pcs.
VACUETTE® Blood Collection Sets + Blood Culture Holder (sterile)				
450186	Blood Collection Set + Blood Culture Holder	tubing length: 19cm 21G x 3/4"	24 pcs.	240 pcs.
Blood Culture Holder (non-sterile)				
450181	Blood Culture Holder		40/bag	320 pcs.

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