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Email : bip@clickTA.com

How to keep a sample "fresh and stored" in the laboratory.

10 rules and some recommendations

1. The procedure is governed by the stability of the constituents of the sample. The most important causes for alterations to the quality of specimen are :

- Metabolism of the blood cells
- Evaporation /sublimation
- Chemical reactions
- Microbiological decomposition
- Osmotic processes
- Effect of light
- Gas diffusion

2. Rapid transport and short storage times improve the reliability of laboratory results.

3. Specimens and samples are preserved longer the cooler they stored (but note exceptions!).

4. Specimens and samples should always be stored in closed vessels (**evaporation!**)

5. The danger of evaporation also exists in refrigerators (condensation of moisture on the cooling elements).


Fig.1 Formation concentration gradients of calcium in plasma sample after rethawing without mixing

6. Storage problems are reduced if disposable sampling systems are used.
7. Separating agents (e.g. gel separators) improve the serum/plasma yields and enable serum to be left in the original tubes above the blood (1).
8. Avoid shaking the sample vessels (pneumatic tube dispatch systems.): risk of hemolysis.
9. Always store sample vessels containing blood vertically; the clotting procedure is accelerated. Label infectious material and handle it with particular care.
10. Label infectious material and handle it with particular care.

8 special rules and some more useful recommendations

1. Avoid storing of whole blood. Information on sensitive analytes is given in the "**list of analytes and leaflets**". Blood samples should reach the laboratory within 45 min of collection in order to ensure that centrifugation and separation of the sample is carried out within 1 hour (2,3,4)
2. Avoid glycolysis to keep glucose, lactate and pH stable. Glycolysis can be avoided by the addition of an inhibitor in conjunction with an anticoagulant (5).
3. Avoid the effect of light otherwise there will be a fall in the values of bilirubin, vitamin C, porphyrins, creatine kinase (CK) and folic acid.
4. Reduce contact with air as far as possible. If this is not done, evaporation/sublimation will result in an apparent increase in the concentration/activity of all non-volatile components. This is particularly the case when the volume of the sample is relatively small and the surface area is relatively large.
5. Whole blood should not be stored in the refrigerator. When urine is cooled, salts may precipitate out of the solution (calcium and magnesium phosphate, uric acid).
6. For certain analytes, the specimens/sample should not be deep frozen. Failure to observe this can result in deviating results for the following analytes:

Table 1 *Examples of blood and urine constituents which should not be stored frozen*

<i>Sample</i>	<i>Analytes</i>
<i>Serum/plasma</i>	<i>Lipoprotein electrophoresis</i> <i>Lipoprotein X</i> <i>Apolipoprotein A-I and B</i> <i>LDL-cholesterol (prevented by the addition of glycerol)</i> <i>Fibrin monomer positive plasma*</i>
<i>EDTA-Blood</i>	<i>Hematology</i>
<i>Urine</i>	<i>IgG</i> <i>Sediment</i> <i>Uric acid (precipitations!)</i>

*Negative test result, prolonged PTT, shortened thrombin time, shortened reptilase time.

7. Correct thawing. A very common source of error is the inadequate mixing of deep-frozen samples after they have thawed. Concentration gradients are produced during thawing as the concentrated solution first melts and then runs down the sides of the vessel (see Fig.1)

After thawing, the sample should therefore be inverted several times, avoiding the formation of foam. Look for undissolved material and, if necessary bring into solution by careful warming.

8. Store sample after analysis in such a way as to permit the confirming of results, checking the identity of sample or performing additional tests for medical or legal reasons:

Table 2 Recommended storage time and conditions for analytical sample

Sample for	Storage time	Temperature
Clinical Chemistry:	1 week	Refrigerator
Immunology:	1 week	Refrigerator
Hematology:	2 day	Room temperature
Coagulation:	1 day	Refrigerator
Toxicology:	6 weeks	Refrigerator
Blood grouping:	1 week (at least)	Refrigerator

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
References

1. Cown and Steels Manual for the Identification of Medical Bacteria Third Edition 1998
2. NCCLS 1990 Quality Assurance for Commercially Prepared Microbiological Culture Media Approved Standard vol 10 No. 4 NCCLS Document M22-A
3. USP XXIV (1999)
4. European Pharmacopoeia III (Addendum 1999)
5. Van Scoy R.E. Culture – Negative Endocar Ohtis–Mayo (lin Proc, 1982, 5713)

Petri dish

- Broad choice of sizes and designs
- Can be filled on semi and fully automated equipment
- Unvented, single and triple vented versions available
- Heavy design available for high temperature agars

Product Number	Product Name
633180	Petri dish, with triple vents 90mmx16 mm
Qty/bag	20pcs/bag
Qty/case	480pcs/Carton

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