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Email : bip@clickTA.com

May I take a coffee (with or without sugar), smoke before blood sampling?

Caffeine

Caffeine is found in many constituents of food ingested daily. Despite its widespread use, the influence of caffeine on various analytes in clinical chemistry has not been investigated in detail. Caffeine inhibits phosphodiesterase and hence cyclic AMP degradation. Cyclic AMP in turn promotes glycogenolysis, there by increasing blood glucose concentrations. In addition, the glucose concentration increases due to gluconeogenesis via epinephrine. Activation of triglyceride lipase leads to a three-fold increase of non-esterified fatty acids¹. Quantification of hormones and drugs bound to albumin is hampered by the fatty acid-induced displacement effect. Three hours after the intake of 250 mg of caffeine, plasma renin activity and catecholamine concentrations have been found to be elevated². Consequently, studies intended to investigate these analytes should take caffeine consumption into account .

Effects of smoking

Smoking leads to a number of acute and chronic changes in analyte concentrations, the chronic changes being rather modest. Smoking increases the plasma/serum

concentrations of fatty acids, epinephrine, free glycerol, aldosterone and cortisol³. These changes occur within one hour of smoking 1-5 cigarettes. Alterations in analytes induced by chronic smoking include blood leukocyte count, lipoproteins, the activities of some enzymes, hormones, vitamins, tumor markers, and heavy metals [Fig 1-1]⁴


Fig.1-1 Deviation of blood analyte concentrations between current smokers and non smokers, chronic effect³

The mechanism underlying these changes has not been fully elucidated. A large number of pyridine compounds, hydrogen cyanide and thiocyanate are found in tobacco smoke. They can account for concentration changes by direct or indirect

effects. Decreased angiotensin converting enzyme activity [ACE] in smokers is believed to result from the destruction of lung endothelial cells with a subsequent reduction in the release of ACE into the pulmonary circulation and/or enzyme inhibition⁵. The extent of changes also depends on the amount, kind of cigarettes/ cigars/pipes and technique of smoking [with or without inhalation]. Moreover, smoking-induced changes are influenced by age and gender⁶.


Fig.1-2 Effect of smoking on different blood analytes caused by smoke constituents^{1,3}

In Fig.1-2 shows the concentrations of cotinine, thiocyanate and carboxy – hemoglobin, used as markers for the qualitative and quantitative assessment of smoking habits. Cotinine has the advantage of having a longer half-life [20-28 hours] than nicotine, the parent compound[12-15 min]⁶

References

- 1.Narayanan S., Physiological variables in blood sampling. Mitt Klin Chem 1993; 24:130-4
- 2.Robertson D., Effect of caffeine on plasma renin activity, catecholamines and blood pressure, New Engl J Med 1978; 298:181-6
- 3.Young DS., Effects of preanalytical variables on clinical laboratory tests. Washington:AACC Press ,1993.
- 4.Haboubi NAA., Serum angiotensin converting enzyme activity in cigarette smokers. Clin Chim Acta 1986; 154: 69-72
5. Startland BE., Selected Pre-analytical Sources of Variation Reference Values in Laboratory Medicine ,1981:127-37
6. Sepkovic DW.Thyroid activity in cigarette smokers, ARCH Intern Med 1984;144:501-3

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VACUETTE® urine system

- ▶ Improved hygiene**
Improved safety for laboratory and hospital personnel through the use of a closed system; contamination is virtually impossible when decanting into a secondary vessel.
- ▶ Safe and convenient transport**
VACUETTE® Urine Tubes are hermetically sealed, sterile, unbreakable and space saving.
- ▶ Pre-analytical standardization**
Errors produced as a result of collection, transportation and storage are reduced to a minimum. Sterile VACUETTE® Urine tubes have an exact dosage vacuum and are guaranteed to meet highest quality standards.
- ▶ Highly efficient and cost saving**
Collection, transportation, centrifugation, strip tests, analysis and storage can all be executed in a single sterile tube.

Instructions for use:

- 1** Insert the urine sample device into the urine sample.
- 2** Push the red rubber seal into the sample device until the needle penetrates the coated cap. The urine will flow automatically into the tube in accordance with the sound defined vacuum.
- 3** Remove the blood test tube from the urine sample device.
- 4** Insert the tube several times after filling. Dispose of the urine sample device and the urine beaker accordingly.

VACUETTE®
one step ahead

The practical solution for standardized pre-analysis

VACUETTE®
Round - base tube (pointed)

VACUETTE®
Control - base tube (pointed)

URINE BEAKER
The hygienic multiplier - sterile, break-proof and leak-proof.

URINE SAMPLE DEVICE
For the convenient sample transfer from the collection vessel.

Recommendations to guarantee a stable sample quality:

- ⌚ In cases where the sample remains in the urine beaker for a longer period of time, the urine sample device should be used to thoroughly mix the sample, thus distributing the sedimentation throughout the sample prior to decanting.
- ⌚ If a sterile urine beaker is used exclusively, after about 1/2 hour bacterial growth occurs, which in turn has an effect on sample quality.
- ⌚ When using urine culture tubes, they should be inverted several times, so as to mix the urine and boric acid homogeneously.
- ⌚ Recommended centrifugal time for pointed-base tubes: 5 Min. / 400 G.

www.vacurette.com

greiner bio one

MiniCollect® Capillary Blood Collection System

- ▶ "Cross Cuts" Cap for filling and collection of the blood sample without the need for opening the tube. Therefore preventing air ingress, hemolysis and possible contamination.
- ▶ Cap with an automatic re-sealing capability
- ▶ Specially designed funnel made of Polypropylene
- ▶ Plastic tube with graduation marks
- ▶ Carrier tube for centrifugation
- ▶ Safety Lancet with various puncture depths
- ▶ Capillary made of specially coated plastic

to collect **to collect** **to collect**

to collect **to collect** **to collect**

to collect **to collect** **to collect**

- ⌚ Flexible rubber cap with "Cross Cuts" which makes de-capping of the tube unnecessary
- ⌚ Maximum safety through a closed system
- ⌚ Suitable for use with a funnel or capillary
- ⌚ Graduation marks to determine exact filling amount
- ⌚ A single tube for blood collection and centrifugation
- ⌚ New improved tube design with free standing ability

MiniCollect®
by VACUETTE

...simple handling with improved safety

Instructions for use

- 1** Clean the puncture area. Remove the lancet protective cover by breaking the cover off.
- 2** Puncture the heel area (venotome) by pressing at the coloured button to perform the puncture.
- 3** Puncture the finger (phlebotomy). Hold the puncture finger between the thumb, pointer and middle fingers. Press in the coloured button to perform the puncture.
- 4** Using slight thumb pressure form a blood drop at the puncture site. With the aid of a capillary or funnel, transfer the blood into the tube. Funnel: Ensure that the funnel handle is positioned over one of the corners of the cap. Capillary: Blood will flow effectively into the tube as a result of the capillary effect.
- 5** After blood collection, gently remove the funnel or capillary from the tube and dispose of. The cap closes automatically. Following blood collection invert the tube several times.

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“ A result is only as good as the specimen. ”

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- ➡ Specimen Collection
- ➡ Transport, Storage of Specimen
- ➡ Specimen Processing and Handling
- ➡ Patient Clinical State

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PHLEBOTOMY TECHNIQUE AND PRACTICES PRE-ANALYTICAL VARIABLES

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4. Levens JK., New millennium, new challenges for coagulation. ADVANCE for Administrator of the Laboratory 1998;7(12):17-20.
5. Jenson R, Fritzma GA., Pre-analytical variables in the coagulation laboratory. DAVANCE for Administrators of the Laboratory 2000;9(7):90-94.
6. National Committee for Clinical Laboratory Standards: Collection, Transport and Processing of Blood Specimens for coagulation testing and performance of Coagulation Assays. Approved Guideline 3rd Ed. December 1998, Document H21-A3, 18, No. 20 p.2-3.
7. Intravenous Nurses Society., Revised Intravenous Nursing Standards of Practice, Standard #33, 1998;21:51-52.
8. Adcock Dm, Kressin Dc. Marlar RA. Effect of 3.2% vs 3.8 % sodium citrate concentration on routine coagulation testing. Am J Clin Pathol 1997;107:150-110.

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