Patient Safety & Specimen Safety that's what we are all about !





SELECT



www.en.kobe-hilo.com

SELECT is an automated collection tube labeler that dispenses correct tube type labeled with barcoded ID as well as texts for patient name and various information vital for tests and clinical managements. SELECT is designed to perform all blood collection sites incuding bedside accessioning, reduce human errors and turn around time.

The best choice for your needs:

- •Holds multiple tube types and sizes in 6 interchangeable stockers.
- Programmable tube stockers.
- •Manual Tube Insert in addition to auto driven stockers.
- •System interfaces with LIS or HIS.
- ·User friendly software and touch screen monitor.
- Customizable bar coded labels including 2D barcodes.
- •Modular design.
- Variety of optional items realize ALL-IN-ONE configuration.

Secondary Printer



Manual Entry for extra tubes

SELECT

User Friendly Software



An Output



Specifications

SELEC

Applicable tubes	12-16mm dia, 75-110mm length Plastic, Rubber, Film or Non capped tubes
Throughput	1.5 sec per a labeling 4 sec. between tubes Only 7 sec. first tube out
Tube types	Up to 6 types in stocker (standard) + manual insert
Number of stocker	6 stockers (standard)
Number of tube/stocker	20 tubes
Label size	50(W) x 30(H)mm (standard)
Labeling position	0-50mm (from top of tube cap)
Operation Panel	LCM Display, 16 alphanumeric x 2 line Keys: STOP/RESTART, FEED, F1, F2
Barcodes	CODE39, CODABAR, EAN128, CODE128, Int 2/5, UPCA/E, EAN 8/13
2D codes	PDF417, QR CODE, MAXI CODE, DATAMATRIX
Fonts	6-30p CG Triumvirate, 16x26 dots US ASCII 8bit, OCR-B, OCR-A, Asian fonts
Print method	Direct Thermal
Print density	8 dots/mm (204 dpi)
Print width	25-54 mm
Print speed	2-4 ips
Interface	RS232C (9600-115200bps), Ethernet (option)
Power	AC100-240V 50/60Hz
power Consumption	0.7A
Weight	37Kg
Dimensions	500 (W) x 490 (D) x 390 (H) mm
Approvals	CE

[OPTIONS]



12" Touch Screen Monitor With Flexible Arm





HILO CO. LTD

* Appearance and specifications are subject to change without prior notice for improvements

Towa Bldg.II #207, 5-8-1 Oishihigashi-machi, Nada-ku, Kobe, Hyogo 657-0043 Japan T: +81 78 414 8305 F: +81 78 414 8306 contact@kobe-hilo.co.jp URL: www.en.kobe-hilo.com